

SINGLE StEP™ 10mL or 4g Reservoir Patented with Fri

Cat# 9452086-100

Size : 100/case

Brand : Thomson Instrument Company

[click here to download a PDF of our SINGLE StEP® Empty Column General information flier](#)

Easy to Use...

At Thomson, we are aware of the need to customize available apparatus to the individual experiments. Our SINGLE StEP® opportunity to fill different sized columns with a variety of sorbents and resins for purification applications. The wide range in are not limited by column size. SINGLE StEP® Empty Columns allow for the simple connection to FPLC/LPLC/MPLC system

1.2L Reservoir for SINGLE StEP® Empty Columns

The Thomson 1.2L Reservoir simplifies gravity filtration protocols simplifying equilibration, sample loading & elution steps. The top of 25g, 40g and 80g SINGLE StEP® Empty Columns allows for easy loading of buffers, solvents and samples minimizing

Image not found or type unknown

Thomson instrument Company is not affiliated with GE or their products.

Key Column Features

- Acceptable for use with Gravity or FPLC/LPLC/MPLC
- Multiple sizes for scales from 10mL-650mL (4g – 300g)
- Durable design for pressures up to 200psi
- Top & bottom connections are standard Luer sizing

Key Reservoir Features

- No mess, screws directly onto a SINGLE StEP® Empty Column
- Easy loading of buffers, solvents and samples