

# SINGLE StEP™ 100mL or 40g Reservoir Patented and F

Cat# 9452092-100

Size : 100/case

Brand : Thomson Instrument Company

You might also be interested by the following products:

Cat#	Description	Size	Price Excl. VAT
VC2571	Vivacell® 250 INSERT; 0.2um; PES Membrane	5 pcs	
VCA250	Vivacell® 250; complete with accessories	1 unit	
9452092-10	SINGLE StEP™ 100mL or 40g Reservoir Patented and Frit	10/case	


[click here to download a PDF of our SINGLE StEP® Empty Column General information flier](#)

## Easy to Use...

At Thomson, we are aware of the need to customize available apparatus to the individual experiments. Our SINGLE StEP® opportunity to fill different sized columns with a variety of sorbents and resins for purification applications. The wide range in are not limited by column size. SINGLE StEP® Empty Columns allow for the simple connection to FPLC/LPLC/MPLC system

## 1.2L Reservoir for SINGLE StEP® Empty Columns

The Thomson 1.2L Reservoir simplifies gravity filtration protocols simplifying equilibration, sample loading & elution steps. The top of 25g, 40g and 80g SINGLE StEP® Empty Columns allows for easy loading of buffers, solvents and samples minimizing



Thomson instrument Company is not affiliated with GE or their products.

## Key Column Features

- Acceptable for use with Gravity or FPLC/LPLC/MPLC
- Multiple sizes for scales from 10mL-650mL (4g – 300g)
- Durable design for pressures up to 200psi
- Top & bottom connections are standard Luer sizing

## Key Reservoir Features

- No mess, screws directly onto a SINGLE StEP® Empty Column
- Easy loading of buffers, solvents and samples