



## SNAP2BHMN0100 | SNAP i.d.® 2.0 Mini Blot Holders (7.5 x 8.4 cm)



The new Mini Blot Holders ensure the membrane is horizontal allowing uniform Antibody distribution. The flat profile allows for compact space allocation for storage, shipping & waste management.

### SNAP2BHMN0100

100/pk

Le prix n'a pas pu être récupéré

La quantité minimale doit être un multiple de

Les prix peuvent être modifiés sans préavis

[Aperçu](#) [Documentation complémentaire](#) [Produits & Applications associés](#)

### Aperçu

[Description](#)  
[Product Information](#)  
[Applications](#)  
[Product Usage Statements](#)  
[Storage and Shipping Information](#)  
[Packaging Information](#)  
[Tableau de caractéristiques principal](#)

### Tableau de caractéristiques principal

| Key Applications                     |                                                                                                                                                                                                                                      |
|--------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| WB                                   |                                                                                                                                                                                                                                      |
| Description                          |                                                                                                                                                                                                                                      |
| Catalogue Number                     | SNAP2BHMN0100                                                                                                                                                                                                                        |
| Trade Name                           | <ul style="list-style-type: none"> <li>Snap i.d.</li> </ul>                                                                                                                                                                          |
| Description                          | SNAP i.d.® 2.0 Mini Blot Holders (7.5 x 8.4 cm)                                                                                                                                                                                      |
| Overview                             | The new Mini Blot Holders ensure the membrane is horizontal allowing uniform Antibody distribution. The flat profile allows for compact space allocation for storage, shipping and waste management.                                 |
| Materials Required but Not Delivered | Vacuum Pump (WP62 or WP61)                                                                                                                                                                                                           |
| Product Information                  |                                                                                                                                                                                                                                      |
| Detection method                     | Chemiluminescence Fluorescence Chromogenic                                                                                                                                                                                           |
| Compatible Blocking Solutions        | <ul style="list-style-type: none"> <li>bløk®-CH Buffer (cat. no. WBAVDCH01)</li> <li>bløk®-FL Buffer (cat. no. WBAVDFL01)</li> <li>bløk®-PO Buffer (cat. no. WBAVDP001)</li> <li>Commercially available blocking reagents</li> </ul> |
| Compatible Detection Reagents        | <ul style="list-style-type: none"> <li>Luminata Western HRP Substrates</li> <li>Commercially available detection reagents</li> </ul>                                                                                                 |
| Compatible Membranes                 | <ul style="list-style-type: none"> <li>PVDF (Immobilon membranes)</li> <li>Nitrocellulose</li> </ul>                                                                                                                                 |
| Applications                         |                                                                                                                                                                                                                                      |
| Application                          | The new Mini Blot Holders ensure the membrane is horizontal allowing uniform Antibody distribution. The flat profile allows for compact space allocation for storage, shipping & waste management.                                   |
| Key Applications                     | <ul style="list-style-type: none"> <li>Western Blotting</li> </ul>                                                                                                                                                                   |

| Applications      |                         |
|-------------------|-------------------------|
| Application Notes | For use on immunoblots. |

  

| Product Usage Statements |                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Usage Statement          | <ul style="list-style-type: none"> <li>Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.</li> </ul> |

  

| Storage and Shipping Information |                  |
|----------------------------------|------------------|
| Storage Temperature              | Room Temperature |

  

| Packaging Information |        |
|-----------------------|--------|
| Material Size         | 100/pk |

[Nous contacter](#)

[Groupe Merck](#) | [Mentions légales](#) | [Conditions d'utilisation](#) | [Respect de la vie privée](#) | [Conditions de vente](#)

© 2018 Merck KGaA, Darmstadt, Allemagne et/ou ses filiales. Tous droits réservés.

© Merck KGaA, Darmstadt, Allemagne, 2014. Toutes les références à Merck désignent Merck KGaA, Darmstadt, Allemagne.