



24833-U Supelco

Capillary Column Nut (fits Agilent injectors and non-MS detectors), hexagonal wrenchtight version

pkg of 2 ea



FDS

Conditionnement - SKUDisponibilité		Prix (EUR)	Quantité
24833-U	Seulement 6 en stock (d'autres en cours d'arrivage) - A PARTIR DE	44.50	0
Commandes Bulk?		AJOUTER AU PANIER	

Propriétés

Related Categories	Analytical/Chromatography , Capillary Column Nut (fits Agilent Injectors and non-MS detectors), hexagonal wrenchtight version , Capillary GC Column Nuts , Column Ferrules, Nuts, and Specialized Hand Tools , GC Accessories , Gas Chromatography (GC) Less...
packaging	pkg of 2 ea

Description

General description
These stainless steel column nuts are replacements for damaged or misplaced original Agilent nuts. Use with *short* ferrules. Hexagonal shape allows tightening with a wrench to keep fingers from getting burned.

Informations Sécurité

RIDADR	NONH for all modes of transport
WGK Germany	3

Documents

Certificat d'Analyse

Entrez le Lot N°

[Devis/ Commande de produits Bulk](#) [FDS](#)

Customers Also Viewed

 22470-U Supelco Ferrule Nut Adapter (fits Agilent injectors and non-MS detectors), hexagonal wrenchtight version pkg of 2 ea	 28034-U Capillary Column Nut (fits Agilent [MS detector] and PerkinElmer®), hexagonal wrenchtight version pkg of 5 ea	 28033-U Capillary Column Nut (fits Varian®) pkg of 1 ea	 22517 Supelco Ferrule Nut Adapter (fits Agilent MS detectors), knurled fingertight version pkg of 2 ea	 22509 Supelco Ferrule Nut Adapter (fits Agilent injectors and non-MS detectors), knurled fingertight version pkg of 2 ea
---	--	--	--	---

Recently Viewed

 24823-U Supelco® M-2A Capillary Ferrule, Long Design I.D. 0.8 mm, 0.53 mm Column I.D., pkg of 10 ea	 24818-U Merlin Microseal™ System Replacement Septum 1 SPME septum	 24816-U Merlin Microseal™ System Replacement Septum 1 General Purpose (3-100 psi) septum	 24817-U Merlin Microseal™ System (fits Varian®) For 1079 injector; 1 nut, 1 inlet adapter, 1 o-ring, and 1 General Purpose (3-100 psi) septum	 24812-U Supelco® M-4 Capillary Ferrule, Short Design pkg of 50 ea, I.D. 0.8 mm, 0.53 mm Column I.D.
---	---	--	--	---

Protocoles et articles

Documentation référencée

[Avez-vous utilisé ce produit dans vos Protocoles et articles? Si oui, cliquez ici.](#)

Produits associés