



XKEM00107 | Pompe à vide manuelle, polypropylène



XKEM00107

1

Le prix n'a pas pu être récupéré

La quantité minimale doit être un multiple de

Les prix peuvent être modifiés sans préavis

[Aperçu](#) [Documentation complémentaire](#) [Produits & Applications associés](#)

Aperçu

Description
Informations physico-chimiques
Informations sur l'emballage

Description	
Référence	XKEM00107
Description	Pompe à vide manuelle, polypropylène
Informations générales	Use polypropylene syringes to dispense liquids into small containers and pressure-filter small volumes through filter holders of filter units with female Luer-Slip inlets. Syringes have male Luer-Slip outlets. An optional two-way vacuum valve, nylon male-to-male Luer adapter and rubber tubing with male Luer inlet converts this standard syringe into a hand vacuum pump. For sterile applications, an autoclavable stainless steel vacuum syringe with stainless steel, two-way valve is available. Applications: Small-Volume Liquid Filtration, Hand Pump for In-Field Filtration with Sterifil Units and Other Filter Holders

Informations physico-chimiques

Pression/Vide	Vide
---------------	------

Informations sur l'emballage

Quantité	1
----------	---

[Nous contacter](#)[Groupe Merck](#) | [Mentions légales](#) | [Conditions d'utilisation](#) | [Respect de la vie privée](#) | [Conditions de vente](#)