

L A N D A N G E R	FICHE TECHNIQUE PRODUIT PRODUCT TECHNICAL SHEET	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-------------------	--	--

Editée le / Publishing date

2022-04-28

Page

1

REFERENCE	DESIGNATIONS - DESCRIPTIONS	CODE EAN
J539104	VALVE MC CULLOCH LARGE 6 DENTS 40X27 MM, REVETEMENT NOIR	3661219308577
	MC CULLOCH WIDE BLADE, 6 TEETH, 40X27MM, BLACK COATING	

Catalogues / catalogs

OrthoVertebre 2019

Chapitre-page

227

DESCRIPTIF DU DISPOSITIF/ PRODUCT DESCRIPTION

L A N D A N G E R

Écarteurs - Retractors - Separadores

Mc Culloch

Valves seules - Blades only - Valvas sólas



Revêtement céramique noir
Ceramic black coating
Revestimiento cerámico negro



30 mm - 1^{1/2}" x 17 mm - 20[#] J539303

40 mm - 1^{3/4}" x 17 mm - 20[#] J539304

50 mm - 2" x 17 mm - 20[#] **J539305**

60 mm - 2^{1/4}" x 17 mm - 20[#] J539306

70 mm - 2^{3/4}" x 17 mm - 20[#] J539307

80 mm - 3^{1/2}" x 17 mm - 20[#] J539308

90 mm - 3^{3/4}" x 17 mm - 20[#] J539309

50 mm - 2" x 20 mm - 45[#] **J539405**

60 mm - 2^{1/2}" x 20 mm - 45[#] J539406

70 mm - 2^{3/4}" x 20 mm - 45[#] J539407




30 mm - 1^{1/2}" x 20 mm - 45[#] J539203

40 mm - 1^{3/4}" x 20 mm - 45[#] J539204

50 mm - 2" x 20 mm - 45[#] **J539205**

60 mm - 2^{1/2}" x 20 mm - 45[#] J539206

70 mm - 2^{3/4}" x 20 mm - 45[#] J539207

80 mm - 3^{1/2}" x 20 mm - 45[#] J539208

90 mm - 3^{3/4}" x 20 mm - 45[#] J539209

100 mm - 4" x 20 mm - 45[#] J539210

120 mm - 4^{3/4}" x 20 mm - 45[#] J539212

30 mm - 1^{1/2}" x 27 mm - 1^{1/10}" J539103

40 mm - 1^{3/4}" x 27 mm - 1^{1/10}" J539104

50 mm - 2" x 27 mm - 1^{1/10}" **J539105**

60 mm - 2^{1/2}" x 27 mm - 1^{1/10}" J539106

70 mm - 2^{3/4}" x 27 mm - 1^{1/10}" J539107

80 mm - 3^{1/2}" x 27 mm - 1^{1/10}" J539108

90 mm - 3^{3/4}" x 27 mm - 1^{1/10}" J539109

100 mm - 4" x 27 mm - 1^{1/10}" J539110

30 mm - 1^{1/2}" x 20 mm - 45[#] J539203

40 mm - 1^{3/4}" x 20 mm - 45[#] J539204

50 mm - 2" x 20 mm - 45[#] **J539205**

60 mm - 2^{1/2}" x 20 mm - 45[#] J539206

70 mm - 2^{3/4}" x 20 mm - 45[#] J539207

80 mm - 3^{1/2}" x 20 mm - 45[#] J539208

90 mm - 3^{3/4}" x 20 mm - 45[#] J539209

100 mm - 4" x 20 mm - 45[#] J539210

120 mm - 4^{3/4}" x 20 mm - 45[#] J539212

30 mm - 1^{1/2}" x 27 mm - 1^{1/10}" J539103

40 mm - 1^{3/4}" x 27 mm - 1^{1/10}" J539104

50 mm - 2" x 27 mm - 1^{1/10}" **J539105**

60 mm - 2^{1/2}" x 27 mm - 1^{1/10}" J539106

70 mm - 2^{3/4}" x 27 mm - 1^{1/10}" J539107

80 mm - 3^{1/2}" x 27 mm - 1^{1/10}" J539108

90 mm - 3^{3/4}" x 27 mm - 1^{1/10}" J539109

100 mm - 4" x 27 mm - 1^{1/10}" J539110

30 mm - 1^{1/2}" x 20 mm - 45[#] J539203

40 mm - 1^{3/4}" x 20 mm - 45[#] J539204

50 mm - 2" x 20 mm - 45[#] **J539205**

60 mm - 2^{1/2}" x 20 mm - 45[#] J539206

70 mm - 2^{3/4}" x 20 mm - 45[#] J539207

80 mm - 3^{1/2}" x 20 mm - 45[#] J539208

90 mm - 3^{3/4}" x 20 mm - 45[#] J539209

100 mm - 4" x 20 mm - 45[#] J539210

120 mm - 4^{3/4}" x 20 mm - 45[#] J539212

30 mm - 1^{1/2}" x 27 mm - 1^{1/10}" J539103

40 mm - 1^{3/4}" x 27 mm - 1^{1/10}" J539104

50 mm - 2" x 27 mm - 1^{1/10}" **J539105**

60 mm - 2^{1/2}" x 27 mm - 1^{1/10}" J539106

70 mm - 2^{3/4}" x 27 mm - 1^{1/10}" J539107

80 mm - 3^{1/2}" x 27 mm - 1^{1/10}" J539108

90 mm - 3^{3/4}" x 27 mm - 1^{1/10}" J539109

100 mm - 4" x 27 mm - 1^{1/10}" J539110

30 mm - 1^{1/2}" x 20 mm - 45[#] J539203

40 mm - 1^{3/4}" x 20 mm - 45[#] J539204

50 mm - 2" x 20 mm - 45[#] **J539205**

60 mm - 2^{1/2}" x 20 mm - 45[#] J539206

70 mm - 2^{3/4}" x 20 mm - 45[#] J539207

80 mm - 3^{1/2}" x 20 mm - 45[#] J539208

90 mm - 3^{3/4}" x 20 mm - 45[#] J539209

100 mm - 4" x 20 mm - 45[#] J539210

120 mm - 4^{3/4}" x 20 mm - 45[#] J539212

30 mm - 1^{1/2}" x 27 mm - 1^{1/10}" J539103

40 mm - 1^{3/4}" x 27 mm - 1^{1/10}" J539104

50 mm - 2" x 27 mm - 1^{1/10}" **J539105**

60 mm - 2^{1/2}" x 27 mm - 1^{1/10}" J539106

70 mm - 2^{3/4}" x 27 mm - 1^{1/10}" J539107

80 mm - 3^{1/2}" x 27 mm - 1^{1/10}" J539108

90 mm - 3^{3/4}" x 27 mm - 1^{1/10}" J539109

100 mm - 4" x 27 mm - 1^{1/10}" J539110

30 mm - 1^{1/2}" x 20 mm - 45[#] J539203

40 mm - 1^{3/4}" x 20 mm - 45[#] J539204

50 mm - 2" x 20 mm - 45[#] **J539205**

60 mm - 2^{1/2}" x 20 mm - 45[#] J539206

70 mm - 2^{3/4}" x 20 mm - 45[#] J539207

80 mm - 3^{1/2}" x 20 mm - 45[#] J539208

90 mm - 3^{3/4}" x 20 mm - 45[#] J539209

100 mm - 4" x 20 mm - 45[#] J539210

120 mm - 4^{3/4}" x 20 mm - 45[#] J539212

30 mm - 1^{1/2}" x 27 mm - 1^{1/10}" J539103

40 mm - 1^{3/4}" x 27 mm - 1^{1/10}" J539104

50 mm - 2" x 27 mm - 1^{1/10}" **J539105**

60 mm - 2^{1/2}" x 27 mm - 1^{1/10}" J539106

70 mm - 2^{3/4}" x 27 mm - 1^{1/10}" J539107

80 mm - 3^{1/2}" x 27 mm - 1^{1/10}" J539108

90 mm - 3^{3/4}" x 27 mm - 1^{1/10}" J539109

100 mm - 4" x 27 mm - 1^{1/10}" J539110

LANDANGER	FICHE TECHNIQUE PRODUIT <i>PRODUCT TECHNICAL SHEET</i>	FAQ 026 INFO Indice : 2 Date : 01/06/2018
-----------	---	---

Éditée le / *Publishing date*

2022-04-28

Page

2

REFERENCE	DESIGNATIONS - <i>DESCRIPTIONS</i>	CODE EAN
J539104	VALVE MC CULLOCH LARGE 6 DENTS 40X27 MM, REVETEMENT NOIR	3661219308577
	MC CULLOCH WIDE BLADE, 6 TEETH, 40X27MM, BLACK COATING	 3 661219 308577

INFORMATIONS REGLEMENTAIRES / *REGULATORY INFORMATION*

Classe CE :

I

Code GMDN :

46774

CE Class :

GMDN code :

PERFORMANCES REVENDIQUEES / *INTENDED PERFORMANCE*

Dispositif médical de bloc opératoire :

REUTILISABLE/REUSABLE

Operating room medical device :

CARACTERISTIQUES SPECIFIQUE / *SPECIFIC FEATURES*

Type matière :

ACIER INOXYDABLE + REVETEMENT / STAINLESS STEEL + COATING

Material type :

Nuance matière :

X20Cr13(1.4021) + AlTiN

Material grade :

Dimensions :

Voir désignations produits

OrthoVertebre 2019

Outer dimensions :

See products description

Chapitre-pages

227

Nettoyage/Pré-désinfection : Standard pour instrumentation de bloc opératoire réutilisable.

Cleaning / Pre-disinfection : *Standard conditions for reusable surgical instruments*

Stérilisation :

La stérilisation doit être effectuée selon un procédé de stérilisation à la vapeur conformément à la norme ISO 17665-1 (procédé sous vide fractionné) avec un stérilisateur conformément à la norme EN 285.

Méthode de stérilisation à la vapeur saturée	Température	Durée	Temps de séchage	Pression absolue
Pré-vide	132 °C (270 °F)	4 minutes	30 minutes	2,9 bars, 2,8 atm, 41,6 psi
Pré-vide	134 °C (273 °F)	3 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	5 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi
Pré-vide	134 °C (273 °F)	18 minutes	35 minutes	3,0 bars, 3,0 atm, 44,1 psi

Sterilization process :

The Sterilization shall be performed according to ISO 17665-1 validated steam sterilisation procedure (fractional vacuum procedure) in an autoclave in compliance with EN 285.

Saturated steam sterilization method	Temperature	Time	Drying time	Absolute pressure
Pre-vacuum	132 °C (270 °F)	4 minutes	30 minutes	2.9 bar, 2.8 atm, 41.6 psi
Pre-vacuum	134 °C (273 °F)	3 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	5 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	18 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi

EMBALLAGE / *PACKAGING*

Conditionnement :

Pochette plastique en polyéthylène (PE).

Packaging :

Polyethylene plastic pouch (PE).