



Product catalog








Hamilton Medical

HAMILTON
MEDICAL

This document reflects your product selection in Hamilton Medical's online product catalog (e-catalog).

Created on: Mar 13, 2023
Country: International

Legend

	Adults		Single use	 The packaging unit for every product is marked within this box icon
	Pediatric		Reusable	
	Neonates		Autoclavable	

We reserve the right to make changes to the online product catalog and prices for the individual articles. The latest and most comprehensive product information is available online.

Accessories



281722

CO2 Mainstream airway adapter

Using the neonatal flow sensor in combination with either a CO2 neonatal single-use (PN 281720) or re-usable (PN 281722) airway adapter, requires the according connector (281803).

- ET ≤ 4 mm
- Dead space: < 1 cc
- Weight: 14.9 grams
- Pressure drop: 0.68 cmH2O @10 LPM

GTIN:	07630002801317
Requires:	281803 Adapter
For ventilator:	HAMILTON-C1/C1 neo/T1/T1 Military/C2/C3/C6/G5/S1





More information:

www.hamilton-medical.com/e-catalog



Manufacturer:

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

+41 (0)58 610 10 20

info@hamilton-medical.com

www.hamilton-medical.com

The products shown here are not available for purchase by the general public. The information provided here is intended for healthcare professionals only. Always read the labels and follow the product's instructions for use. Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. All images are used for illustrative purposes only and may not accurately represent the product or its use. For all proprietary trademarks (®) and third-party trademarks used by Hamilton Medical AG see www.hamilton-medical.com/trademarks. ©2022 Hamilton Medical AG. All rights reserved.