

## Spike SARS CoV-2 DELTA mutant mRNA (mRNA encoding spike protein of SARS-CoV-2 DELTA mutant)

### Description

Ready-to-use stabilized Spike SARS CoV-2 DELTA mutant (L452R; T478K; D614G; P681R; D950N) mRNA.


Concentration: 1.0mg/mL in 1mM Sodium Citrate, pH6.4.

mRNA length: 4036 nt. MW: 1300873 g/mol.

Spike SARS CoV-2 DELTA mRNAs have been designed to produce high expression level of DELTA mutant spike protein of SARS-CoV-2 virus. OZB mRNAs are produced by *in vitro* transcription. mRNAs are stabilized at the 5' end by modified nucleotides capping (Cap1) and contain a poly(A) tail at the 3' end. Sequences have been optimized to yield improved stability and performance. DELTA Spike protein of SARS-CoV-2 virus mRNA (ref# **MRNA36**) does not bear any additional nucleotide modifications while (ref# **MRNA37**) is modified with 5-methoxyuridine (5moU) to reduce innate immune responses.


(ref# MRNA36):

*Mature mRNA (unmodified nucleotides) with cap1 and polyA tail*



(ref# MRNA37):

*Mature mRNA (fully modified moU) with cap1 and polyA tail*



### Kit contents

**Spike SARS CoV-2 DELTA mRNAs-20:** 20 µg of mRNA unmodified or moU

**Spike SARS CoV-2 DELTA mRNAs-100:** 100 µg of mRNA unmodified or moU

**Spike SARS CoV-2 DELTA mRNAs-1000:** 1 mg of mRNA unmodified or moU

### Storage

Spike SARS CoV-2 DELTA mutant mRNAs must be stored at -80°C. We recommend to aliquot the mRNA solution for a better storage.

### Applications

For the COVID-19 vaccines, scientists developed synthetic mRNA in a lab that instructs cells to produce the distinctive spike protein from the SARS-CoV-2 virus. The immune system then targets and destroys these

foreign spike proteins. If the body encounters the real virus at a later time, the body's immune system will already be prepared to fend it off again. This mRNA encodes for the DELTA mutant Spike protein of the SARS-CoV-2 virus and bears five mutations. Spike protein is used as antigen for immunization and biochemical studies. **For Research Use Only. Not for use in humans. Not for use in diagnostic or therapeutic purposes.** Spike mRNAs resemble fully matured mRNAs with 5'cap1 structure and 3' polyA tail, therefore ready to be translated by the ribosome. mRNA transfection provides several advantages over plasmid DNA (pDNA) delivery. It does not require nuclear uptake for being expressed since translation of mRNA occurs directly into cytoplasm. Indeed, nuclear delivery (transport through nuclear membrane) is one the principal barriers for transfecting slow or non-dividing cells and consequently, mRNA transfection is particularly attractive for such purpose. This approach presents also the advantage of being non-integrative which is particularly appealing for stem cells, regenerative medicine or vaccine fields. Contrary to pDNA, mRNA cannot lead to genetic insertion causing mutations. Moreover, the protein expression from the mRNA is promoter-independent and faster than with DNA. For transfection we recommend RmesFect™ (#RM21000) and RmesFect™ Stem (#RS31000).

### Related Products

Ref	Description
RM21000	RmesFect™ transfection reagent 1mL
RS31000	RmesFect™ Stem transfection reagent 1mL
MRNA41/42	mRNA OVA (unmodified or moU)
MRNA34/35	mRNA spike SARS Cov2 (E484K; N501Y)

Discover the complete list of mRNA at: [www.ozbiosciences.com](http://www.ozbiosciences.com)  
**Custom mRNAs are also available now!**

### Contact Us

Feel free to contact us for all complementary information and remember to visit our website to stay informed on the latest breakthrough technologies and updated on our complete product list. ([www.ozbiosciences.com](http://www.ozbiosciences.com)). For bulk, please contact us: [order@ozbiosciences.com](mailto:order@ozbiosciences.com)

Technical questions: [tech@ozbiosciences.com](mailto:tech@ozbiosciences.com)  
Order: [order@ozbiosciences.com](mailto:order@ozbiosciences.com)  
General Information: [contact@ozbiosciences.com](mailto:contact@ozbiosciences.com)

#### OZ Biosciences SAS

163 avenue de Luminy  
Case 922, zone entreprise  
13288 Marseille cedex 09 - FRANCE  
Ph: +33 (0) 486 948 516  
Fax: +33 (0) 486 948 515  
[contact@ozbiosciences.com](mailto:contact@ozbiosciences.com)  
[order@ozbiosciences.com](mailto:order@ozbiosciences.com)

#### OZ Biosciences INC

4901 Morena Blvd,  
Suite 901  
San Diego CA 92117 - USA  
Ph : + 1-858-246-7840  
Fax : + 1-855-631-0626  
[contactUSA@ozbiosciences.com](mailto:contactUSA@ozbiosciences.com)  
[orderUSA@ozbiosciences.com](mailto:orderUSA@ozbiosciences.com)