

SUB Aqua Pro - 2L, 2L shallow, 5L, 12L, 18L, 26L, 34L and dual

NEW

- | | |
|--------------------|---|
| PERFORMANCE | › Set and Forget™ technology - fast heat-up, accurate temperature control |
| PRECISION | › Excellent temperature stability - $\pm 0.2^{\circ}\text{C}$ |
| QUALITY | › 3 year warranty. UK design and manufacture |
| PROTECT | › User-defined high-temperature alarm & cut-out - protects samples |
| CONVENIENT | › Drain tap, tray and lid included |
| OPTIMISED | › Heater mat technology - optimised work space and easy to clean |
| VALUE | › Advanced dry start and run dry protection - prevents costly service repairs |
| ENHANCED | › 3 programmable presets, countdown timer, 2 point calibration, |



3
year warranty



Applications

- University - research, sample preparation
- Healthcare - warming samples
- Industry - QA/QC samples
- Pharma/biotech - sample preparation

**Ideal choice for University research
and Pharma/Biotech**

Specifications



	SAP2	SAP2S	SAP5	SAP12	SAP18	SAP26	SAP34	SAPD
								
	2.5 kg l: 200mm w: 185 mm h: 200 mm	3 kg l: 215 mm w: 335 mm h: 150mm	3 kg l: 215 mm w: 335 mm h: 200 mm	6 kg l: 380 mm w: 360 mm h: 225 mm	9 kg l: 590 mm w: 335 mm h: 275 mm	9 kg l: 590 mm w: 335 mm h: 275 mm	14.5 kg l: 770mm w: 335 mm h: 370 mm	9 kg l: 380 mm w: 545 mm h: 225 mm
Capacity	2L	2L shallow	5L	12L	18L	26L	34L	5L & 12L
Temperature range	ambient +5 to 99°C							
Temperature display and setting resolution	0.1°C							
Stability (DIN 12876-3 @70 °C)	± 0.2°C							
Working dimensions (mm)	125 x 140 x 115	145 x 290 x 300	145 x 290 x 115	315 x 290 x 115	495 x 290 x 115	495 x 290 x 165	630 x 290 x 160	145 x 290 x 115 & 315 x 290 x 115
User defined calibration points	2							
Drain tap	-	-	-	●	●	●	●	●
Heater power 230V/120V (kW)	0.25/0.25	0.35/0.35	0.35/0.35	0.8/0.8	1.4/1.05	1.4/1.05	1.8/1.3	1.15/1.15
Supply voltage (V)	120 or 230 (50-60Hz)							

Essential accessories (full range available on website)


Flat lid with ring sets

Variable hole diameter to accommodate tall vessels whilst reducing evaporation



Polypropylene spheres

Useful alternative to a lid, minimises evaporation and heat loss whilst allowing easy access to vessels in the bath; particularly useful for tall vessels



Test tube racks

Stainless steel - range of sizes available for 0.5 to 1.5 ml micro-tubes or 10 to 30 mm tubes



Raised shelves

Stainless steel - shelf covers half area of bath. Reversible, allows two shelf depths



Alternative water bath ranges:



JB Nova - general purpose

- Includes clear polycarbonate lid
- 4 bath sizes: 5L, 12L, 18L, 26L
- Drain tap on 12L, 18L and 26L baths

JB Academy - basic

- Simple and practical, great value
- 3 bath sizes: 5L, 12L, 18L



Call us:

+ 44 (0) 1763 260 811

Email:

sales@grantinstruments.com