

FICHE TECHNIQUE

Référence	255008
Désignation	Plasmid maxiPREP Kit



Description :

Sample: Up to 200 ml bacterial cells
Yield: Up to 850 µg of plasmid
Endotoxin Value: <0.003 EU/ug

Features

Safe: Eliminates the use of phenol, chloroform, ethidium bromide, and cesium chloride, thus minimizing the exposure to and disposal of hazardous materials.
Time saving: Complete the process in less than 40 minutes.

Unité de vente :

Size: 10 rxns
