


"Bio-Adembeads Antibodies : Very easy selection"



Indirect IgG/IgM Selection : High flexibility

- **Rapid procedure:** capture of any type of cells from different species in less than 5 min
- **High specificity:**
 - all our coated anti-species antibodies are Fc specific
 - no cross reaction
- **Simple:**
 - all frequently used antibodies can be recognized by Bio-Adembeads Antibodies
 - Compatible with flow cytometry and microscopic analysis



Ordering Information : Bio-Adembeads Antibodies

Product	Size	Code Number
Goat Anti-Mouse IgG	2-10-50 ml	G05030-G05031-G05032
Goat Anti-Rat IgG	2-10-50 ml	G05130-G05131-G05132
Goat Anti-Rabbit IgG	2-10-50 ml	G05230-G05231-G05232
Goat Anti-Human IgG	2-10-50 ml	G05330-G05331-G05332
Mouse Anti-biotin	2-10-50 ml	G05430-G05431-G05432
Goat Anti-Mouse IgM	2-10 ml	M05030-M05031
Goat Anti-Rat IgM	2-10 ml	M05130-M05131
Goat Anti-Rabbit IgM	2-10 ml	M05230-M05231
Goat Anti-Human IgM	2-10 ml	M05330-M05331

Diameter
300 nm

ademtech



ISO 9001 Certified

ADEMTECH SA

Parc BIOGALIEN

27 Allée Charles Darwin - 33600 Pessac - France

Tel: +33 (0)5 57 02 02 01 - Fax: +33 (0)5 57 02 02 06

E-mail: ademtech@ademtech.com - www.ademtech.com