

TAE 25X
Tris Acetate EDTA buffer solution
25X concentrated
Ref. EU0200

Suitable for molecular biology, immunology and nucleic acid methodologies

Buffer composition: - Tris(hydroxymethyl)aminomethane: 1 M
(25X concentrated) - Acetate (acetic acid): 0.5 M
- EDTA: 25 mM

Physical state: Ultra-pure reagents dissolved in ultra-pure deionized water (18.2 MΩ).
Liquid filtered 0.2 µm.
Stabilizer: none
Preservative: none
This product is a 25X concentrated stock solution and should be diluted with deionized water to 1X.

pH: 8.0-9.0 (20°C)
Conductivity: 15-20 mS/cm (20°C)
Resistivity: 60-70 Ωcm (20°C)
Sal: 9-10 (20°C)

Storage: Store at room temperature
Protect from light
 Carefully lock containers and **upright** store, in order to prevent any running out.
Refermer soigneusement les récipients et stocker **debout**, pour éviter tout écoulement.

Handling/safety: Avoid contact with skin, eyes and clothing.



Attention

H315 Provoque une irritation cutanée.

H319 Provoque une sévère irritation des yeux.

P305+P351+P338 EN CAS DE CONTACT AVEC LES YEUX: rincer avec précaution à l'eau pendant plusieurs minutes. Enlever les lentilles de contact si la victime en porte et si elles peuvent être facilement enlevées. Continuer à rincer.