

TRIS-EDTA solution
100X concentrated
Ref. EU0960

Size:	100 ml 500 ml 1 l
Form:	Tris amino ultra pure: 1M EDTA disodium salt: 100mM Prepared with ultra-pure water (18.2 MΩ). Filtered using validated 0.2µm membrane filters
pH_(100X):	8.0 _(HCl)
Preservative:	none
Dilution:	This product is a 100X concentrated stock solution and should be diluted with distilled water before use.
Storage:	Store at room temperature.  Carefully lock containers and upright store, in order to prevent any running out. Refermer soigneusement les récipients et stocker debout , pour éviter tout écoulement.
Handling:	Avoid contact with skin, eyes and clothing. After eye contact: rinse opened eye for several minutes under running water.