









# ImmPRESS®-AP Horse Anti-Rabbit IgG Polymer Detection Kit, Alkaline Phosphatase

## MP-5401

[Product Images](#)





## Short Description

---

The ImmPRESS polymerized reporter enzyme staining system uses novel conjugation and micropolymer chemistries to create a highly sensitive, ready-to-use, one-step detection system. This unique micropolymer of highly active enzyme is attached to our affinity purified secondary antibodies, producing reagents with outstanding sensitivity and low background.

### Features:

- Enhanced accessibility to nuclear and membrane antigens
- Ready-to-use in a convenient dropper bottle
- Shorter assay time
- Simplified multiple labeling

### Kit Contents:

- ImmPRESS-AP (alkaline phosphatase) Polymer Anti-Rabbit IgG Reagent (ready-to-use)
- 2.5% Normal Horse Serum for blocking (ready-to-use)

The ImmPRESS (15 ml) Kit will stain approximately 75-150 sections and the ImmPRESS (50 ml) Kit will stain approximately 250-500 sections.

## Additional Information

---

|                     |                                                                           |
|---------------------|---------------------------------------------------------------------------|
| Unit Size           | 15 ml, 50 ml                                                              |
| Applications        | Immunohistochemistry / Immunocytochemistry, In situ hybridization, ELISAs |
| Technology          | Micropolymer Reagents                                                     |
| Detection Enzyme(s) | Alkaline Phosphatase                                                      |
| Conjugate           | Micropolymer AP                                                           |
| Target Species      | Rabbit                                                                    |
| Format              | Ready-to-Use                                                              |

### ImmPRESS®-AP Horse Anti-Rabbit IgG Polymer Detection Kit, Alkaline Phosphatase



#### Features:

- High sensitivity
- Low background
- Non-biotin
- One-step detection
- Micropolymer technology
- Enhanced accessibility to nuclear and membrane antigens
- Ready-to-use in a convenient dropper bottle
- Shorter assay time
- Simplified multiple labeling

The ImmPRESS (15 ml) Kit will stain approximately 75-150 sections.

### ImmPRESS®-AP Horse Anti-Rabbit IgG Polymer Detection Kit, Alkaline Phosphatase



#### Features:

- High sensitivity
- Low background
- Non-biotin
- One-step detection
- Micropolymer technology
- Enhanced accessibility to nuclear and membrane antigens
- Ready-to-use in a convenient dropper bottle
- Shorter assay time
- Simplified multiple labeling

The ImmPRESS (50 ml) Kit will stain approximately 250-500 sections.



